


Determination of the Ideal Gas Constant (R)

In this lab we are going to experimentally measure the ideal gas constant (R of $PV = nRT$). To do this we will be utilizing the reaction of magnesium with hydrochloric acid. The equation for the reaction is as follows.



Equipment needed

600 mL beaker

rubber stopper with hole

10 mL graduated cylinder

Temperature probe (computer to run it)

1. Obtain a small piece of magnesium ribbon and measure its length in **centimeters** as accurately as you can (to the hundredths place). Record this length in the data table under trial 1.
2. **Do not push the magnesium ribbon all the way into the stopper.** Wedge part of the ribbon into the skinnier end of the rubber stopper, just enough so it won't easily fall through.
3. Fill your 600 mL beaker $\frac{3}{4}$ of the way with tap water.
4. Fill the graduated cylinder all the way up with 1 M HCl (so that it is overflowing).
5. This step must be done quickly!! Place the rubber stopper on the graduated cylinder. Place one finger over the hole in the stopper and flip the graduated cylinder upside-down into the 600 mL beaker. You will see a reaction occurring and gas forming. It is fine if the magnesium floats up into the graduated cylinder.
6. Place your temperature probe in the beaker near the mouth of the graduated cylinder.
7. Wait until the reaction completes (until the bubbling stops).
8. Convert the temperature to Kelvin and record that in your data table.
9. Tap the side of the cylinder a few times to dislodge and air bubbles that may be sticking.
10. Hold the graduated cylinder so that the water level inside the graduated cylinder and the water level of the beaker are the same (this will make the pressure inside the graduated cylinder the same as the pressure in the room). Accurately, read the volume of the gas in **milliliters** by reading the volume in the graduated cylinder (to the tenths place for these graduated cylinders).
11. Convert this volume to Liters and record it in your data table
12. Record the pressure in the room in kilopascals.
13. Dispose of the excess HCl in the waste container under the hood.
14. Repeat this lab two more times and record your data in the column trial 2 and trial 3 respectively.

Calculations

1. The mass of 100.0 cm of the magnesium ribbon used will be written on the board for you. Use this mass as a conversion factor convert your length of magnesium ribbon into grams of magnesium ribbon.
2. Convert your mass of magnesium to moles of magnesium ribbon, using the atomic mass of magnesium. This is the same as the number of moles of H₂ gas.
3. For each trial calculate a value for R by plugging all experimental values into the Ideal Gas Law equation. (check your units, volume must be in liters and temperature in Kelvin)
4. Take the average of all three values of R.
5. Calculate the percent error for your observed value (your average) against the accepted value of R (8.31 L kPa/mol K)

Mass of 100.0 cm of magnesium ribbon _____

	Trial 1	Trial 2	Trial 3
Length magnesium ribbon			
Mass of magnesium used			
Moles of Mg/H ₂ gas (same number)			
Volume of H ₂ gas			
Pressure of room			
Temperature			
Value of R			

(show work for calculation of R)

Average Value of R

Calculate the percent error for the value of R

$$\left| \frac{(\text{observed} - \text{actual})}{\text{actual}} \right| \times 100$$

